

TAPING SYSTEMS CATALOG



**PowerFlex® + PowerTape®  
VictoryTape™**  
THE ULTIMATE TAPING SYSTEM





# PowerFlex® + PowerTape® VictoryTape™

## THE ULTIMATE TAPING SYSTEM

### Our Philosophy

*Taping and wrapping is a skilled art form. Every taping procedure applied by a skilled professional should be a masterpiece...a work of art that an athlete can trust.*

Andover strongly feels our cohesive taping systems are a major performance advancement when compared to the age-old foam pre-wrap and cloth tape methods. The no-stretch, water-resistant, all-day performance you get from PowerFlex® combined with PowerTape® or VictoryTape™ delivers unparalleled support in a taping system.


New products and new methods take time to embrace, it requires an open mind and willingness to break with tradition. Some trainers have still not tried the Andover Ultimate Taping System approach. The Andover team would be happy to discuss samples or any questions you may have. College or professional teams are eligible for a custom presentation by Ron O'Neil. For more information please call 800.432.6686 or visit [andoverhealthcare.com](http://andoverhealthcare.com)

### Our History

Andover Healthcare is the leading manufacturer of cohesive bandages and tapes for almost 40 years and currently supplies the sports medicine, healthcare, and animal health industries with high quality products made in the USA. The company was founded in Massachusetts in 1976, and now distributes its patented products around the world.

### Timeline

- 1995: Andover launches PowerFlex®.
- 2005: Andover launches PowerTape® and the Ultimate Taping System is created.
- 2007: Andover launches PowerFast™ & PowerFlex® AFD.
- 2008: Andover launches PowerSpeed™.
- 2012: The Ultimate Taping System gets more economical with the launch of VictoryTape™.



**Taping Methods Video Textbook** is an educational DVD featuring unique taping methods by Ron O'Neil, ATC. Includes preventative and post-injury taping instruction and taping techniques for ankle and lower leg, knee, upper extremity and more. Great for use as a class supplement.

# PowerFlex®

## SELF-ADHERING Sports Wrap

PowerFlex® is the leader in cohesive sports wraps for its superior performance, wide variety of applications and also for its compatibility as the foundation of Andover's Ultimate Taping System. PowerFlex® is a patented, synthetic, breathable cohesive bandage that can be used alone or as the base layer with cohesive PowerTape®/VictoryTape™ to support muscles and joints. PowerFlex® is sweat and water-resistant and will not slip, slide or migrate for all-day performance. PowerFlex® sticks to itself, not skin or hair, to provide 23 lbs. of tensile strength, the most of any sports wrap on the market. Tears by hand and easy to apply.

- #1 Wrap in Pro Football & College Sports
- Sweat & water-resistant
- All-day performance, won't slip, slide or migrate
- Eliminates the need for pre-wrap and adhesive spray
- Helps prevent injuries
- Tears by hand - patented technology
- Made in the USA



Catalog #	Size	Rolls/Case
3710XX-048	1" x 6 yds.	48
3715XX-032	1.5" x 6 yds.	32
3720XX-024	2" x 6 yds.	24
3725XX-016	2.75" x 6 yds.	16
4725XX-048	2.75" x 6 yds.	48
4730XX-048	3" x 6 yds.	48
3730XX-016	3" x 6 yds.	16
3740XX-012	4" x 6 yds.	12
3760XX-008	6" x 6 yds.	8

**Colors Available:** Red (RD), Blue (BL), Green (GR), Yellow (YL), White (WH), Black (BK), Purple (PU), Neon Green (NG), Neon Pink (NP), Light Blue (LB), Navy (NV), Tan (TN), Teal (TE), Maroon (MR), Orange (OR) and Grey (GY).

*PowerFlex® also available non-latex (White).*

**Team Color Prints Available:** Black Stripe on Orange (TO), Black Stripe on Red (TR), Green Stripe on Yellow (TY), Blue Stripe on Grey (TG), Black Stripe on Neon Pink (TP), Blue Stripe on White (WB), Purple Stripe on White (WP), Green Stripe on White (WG), Orange Stripe on White (WO), Red Stripe on White (WR).



U.S. PATENT NO. 5,762,623  
EUROPEAN PATENT NO. 0839222

# PowerTape<sup>®</sup>

## SELF-ADHERING Sports Tape

PowerTape<sup>®</sup> is the no-stretch, water-resistant cohesive bandage that when combined with PowerFlex<sup>®</sup> creates the Ultimate Taping System. PowerTape<sup>®</sup> is a cohesive athletic tape that sticks only to itself or to PowerFlex<sup>®</sup> creating the strongest combination taping system in existence. PowerTape<sup>®</sup> contains less than 1% stretch - eliminating the “give or slip” that exists when using traditional athletic tape. PowerTape<sup>®</sup> and PowerFlex<sup>®</sup> mold together, creating a taping system with better conformability and comfort. The PowerTape<sup>®</sup>/PowerFlex<sup>®</sup> Taping System is sweat and water-resistant and remains supportive and stable throughout games and practices, even in adverse playing conditions.

- 41 lbs. of tensile strength
- Replaces traditional cloth athletic tape
- Hand-tear cohesive tape that sticks to itself or PowerFlex<sup>®</sup>, not skin or hair
- Less than 1% stretch - eliminates the “give or slip” that is found with cloth tape
- Sweat & water-resistant
- All-day performance
- Excellent conformability and “lay in”
- Made in the USA

Catalog #	Size	Rolls/Case
ACP130XX-1”	1” x 15 yds.	48
ACP130XX-1.5”	1.5” x 15 yds.	32
ACP130XX-2”	2” x 15 yds.	24

**Colors Available:** Tan (TN), Red (RD), Blue (BL), Green (GR), White (WH), Black (BK), Navy (NV), Orange (OR).

*Also available non-latex - must use with non-latex PowerFlex (white only).*

When placing an order, please state catalog # and color. Example ACP130BL. XX = Color Code.



**STEP 1:**  
Apply PowerFlex directly to the skin with proper un-wind tension.



**STEP 2:**  
Cover PowerFlex with PowerTape to provide the strongest taping system in existence.

# VictoryTape™

## ECONOMICAL + EASY TO HAND-TEAR

VictoryTape™ is the economical, easy to hand-tear, self-adhering tape that still delivers 90% of the performance of PowerTape® but at a lower cost. Apply over PowerFlex® for a strong, long-lasting taping system. VictoryTape™ contains less than 1% stretch - eliminating the “give or slip” that exists when using traditional athletic tape. VictoryTape™ and PowerFlex® mold together creating a taping system with better conformability and comfort. The system is sweat and water-resistant, remaining supportive and stable throughout games and practices, as well as in adverse playing conditions.

- **38 lbs. of tensile strength**
- **Replaces traditional cloth athletic tape**
- **Easy hand-tear cohesive tape that sticks to itself or PowerFlex®, not skin or hair**
- **Less than 1% stretch – eliminates the “give or slip” that is found with cloth tape**
- **Sweat and water-resistant**
- **All-day performance**
- **Excellent conformability and “lay in”**
- **Made in the USA**

Catalog #	Size	Rolls/Case
ACP120XX-1”	1” x 15 yds.	48
ACP120XX-1.5”	1.5” x 15 yds.	32
ACP120XX-2”	2” x 15 yds.	24

**Colors Available:** Tan (TN), Red (RD), Blue (BL), Green (GR), White (WH), Black (BK), Navy (NV), Orange (OR).

When placing an order, please state catalog # and color. Example ACP120RD. XX = Color Code.

### What is the difference between PowerTape® & VictoryTape™?

- PowerTape® is the original and still gold standard of cohesive tapes with 41 lbs. of tensile strength and less than 1% stretch.
- VictoryTape™ is a more economical, easy to hand-tear version of PowerTape®. You get 38 lbs. of tensile strength with the same no-stretch, water-resistant performance of PowerTape®.



PATENT PENDING

# PowerSpeed™

STRONGER THAN POWERFLEX

## THE ULTIMATE SPEED TAPE

PowerSpeed™, a cohesive athletic wrap designed for speed taping or for a stronger spat delivers maximum support for a modified taping system. It won't slip/migrate from sweat or water, and tears easily by hand for quick application.

- Designed and developed for speed taping
- Maximum support in a modified taping system
- Strongest cohesive strapping product available
- Sweat & water-resistant for all-day performance - won't slip, slide or migrate
- Eliminates the need for pre-wrap and adhesive spray
- Tears by hand - patented technology
- Made in the USA



**Uses:**

Speed taping, spattng, ankle injury prevention, supports sprains and strains, provides injury prevention and post-injury stability for inversions, eversions and high ankle sprains.

Catalog #	Size	Rolls/Case
PS3720XX-024	2" x 6 yds.	24
PS3725XX-016	2.75" x 6 yds.	16
PS3730XX-016	3" x 6 yds.	16

Colors Available: White (WH), Black (BK).

U.S. PATENT NO. 5.762.623  
EUROPEAN PATENT NO. 0839222

# PowerFlex AFD

## ABSORBENT FOAM DRESSING

PowerFlex® AFD is an all-in-one foam pad dressing and cohesive bandage that delivers instant application mid-game. A highly absorbent foam pad is attached to strong, flexible PowerFlex® for quick and easy application. A waterproof film layer between the PowerFlex® and foam pad prevents leakage and contamination. Fast, effective treatment for wounds, lacerations and turf burns. Available sterile.

Catalog #	Size	Rolls/Case
3722XX-016	2" x 2.5 yds.	16
3744XX-008	4" x 2.5 yds.	8
3722(S)	2" x 2.5 yds.	16
3744(S)	4" x 2.5 yds.	8



Colors Available: White (WH), Black (BK). Sterile (S) = Sterile Product/White only.

U.S. PATENT NO. 5.762.623 PATENT PENDING  
EUROPEAN PATENT NO. 0839222

# PowerFAST™

## FLEXIBLE ADHESIVE STRETCH TAPE

PowerFAST™ is a strong, flexible pressure sensitive adhesive stretch tape made of synthetic fabric. Thin and comfortable, PowerFAST™ won't absorb water and retains elasticity to provide long lasting support.

- Thin and comfortable – no bulk
- Great adhesion yet won't pull hair
- Made in the USA

### Uses:

Ankle injury prevention, supports sprains and strains, provides injury prevention and post-injury stability.



Catalog #	Size	Rolls/Case
4210XX-048	1" x 5 yds.	48
4215XX-032	1.5" x 5 yds.	32
4220XX-024	2" x 5 yds.	24
4230XX-016	3" x 5 yds.	16
4310XX-048	1" x 7.5 yds.	48
4315XX-032	1.5" x 7.5 yds.	32
4320XX-024	2" x 7.5 yds.	24
4330XX-016	3" x 7.5 yds.	16

**Colors Available:** White (WH), Black (BK), Red (RD), Blue (BL), Green (GR), Grey (GY).


When placing an order, please state catalog # and color. Example 4220WH. XX = Color Code. U.S. PATENT NO. 5.762.623 EUROPEAN PATENT NO. 0839222

# Moleskin

## High Quality Protection

Flannel with a napped cotton backing. Use it to prevent blisters, corns, calluses and chafing in athletic footwear.

Catalog #	Size	Rolls/Case
<b>#40 Extra Heavy Moleskin</b>		
040-12-012	12" x 5 yds.	12
040-09-012	9" x 5 yds.	12
040-03-016	3" x 25 yds.	16
040-02-024	2" x 25 yds.	24
<b>#41 Heavy Moleskin</b>		
041-12-012	12" x 5 yds.	12
041-09-012	9" x 5 yds.	12
041-03-016	3" x 25 yds.	16
041-02-024	2" x 25 yds.	24





**PROUDLY MADE IN THE USA**



ANDOVER HEALTHCARE, INC.  
9 FANARAS DRIVE, SALISBURY, MA 01952  
PHONE: 978.465.0044 TOLL FREE: 800.432.6686  
FAX: 978.462.0003  
[www.andoverhealthcare.com](http://www.andoverhealthcare.com)