

**NEW!**

# CoFlex<sup>®</sup> + UMAFD

## UNIVERSAL MEDIC ABSORBENT FOAM DRESSING

### All-in-One Trauma Bandage System

#### Features:

- High absorption and retention of fluids
- All-in-one dressing/cohesive bandage
- Additional foam/packing in core
- Waterproof film layer
- Non-adherent wound pad
- Flexible compression
- Water resistant cohesive bandage
- Expands to conform to wound bed
- Easy-open, vacuum-sealed packaging
- Able to withstand high temperatures
- Up to 140lbs. tensile strength




#### Benefits:

- Quickly secures wound
- Self adhering - holds without clips or fasteners
- Each pad absorbs up to 133 cc's of fluid (8.2 total oz.)
- Prevents leakage and contamination
- Protects wound area, easy to remove
- Helps stop bleeding, conforms to body parts
- Stays in place despite adverse conditions
- Helps stop bleeding and reduces leakage
- Small and easy to pack; rolls won't crush
- Heat tested at 130°F for 90 days
- Provides strong, direct pressure without breaking

### CoFlex UMAFD Construction

U.S. Patent Nos. 5,762,623 6,156,424 Patent Pending  
European Patent No. 0839222



NSN #: 6510-01-571-9729

DAPA #: SP0200-08-H0064

FSS #: V797P-4908A

#### ORDERING INFORMATION

Catalog#	Size	Rolls/Case
9430UM-01	4" x 3 yds.	1
9430UM-018	4" x 3 yds.	18
9430UM-100	4" x 3 yds.	100

Proudly Made in USA



Andover Healthcare, Inc.

9 Fanaras Drive, Salisbury, MA 01952 Toll Free: 800-432-6686 Phone: 978-465-0044 Fax: 978-462-0003 [www.andoverhealthcare.com](http://www.andoverhealthcare.com)